

# PlanTAG

---

Full arch cases • A dual scan protocol

# 1 Final teeth positions for the final restoration must be confirmed.

1.1 If the patients existing denture replicates this it may be used. Make relines prior to the CT scan.



1.2 If not , a diagnostic denture should be made.


2

## 6 radiopaque markers are placed: 4 into the buccal, 2 into the lingual/palatal.

The holes are filed with a radiopaque material, x-ray markers or (e.g.gutta percha, or barium sulfate mixed with acrylic 50/50).The markers should be spread around the arch .T he markers may be placed in the buccal or lingual phalange.




Dual scan  
appliance  
with markers  
attached



We are recommended to share to a CT center.

3 A CT scan is taken first of the denture by itself at 2-3 Ma.

It should be reviewed to sure all radiopaque markers and the denture can be seen.



DENTURE  
INSIDE  
MOUTH!!!

RADIATION  
LEVEL 9-12 Ma

4 A CT scan with the patient wearing the denture is now taken at the conventional settings of 9-12 Ma.

The maxilla and mandible should be separated.



DENTURE  
INSIDE  
MOUTH!!!

RADIATION  
LEVEL 9-12 Ma

tag dental



**T.A.G. Medical Products Corporation Ltd**  
Kibbutz Gaaton 2513000 Israel Tel: +972 (0)4 9858400  
[www.tag-med.com](http://www.tag-med.com)

---

T.A.G. Medical Products Corporation Ltd. is ISO 13485  
and CAN/CSA ISO 13485 certified by LNE

T.A.G. Medical Products Corporation Ltd. Products are  
CE marked in accordance with Directive 93/42/EEC