



COMPACT MASSIF SURGICAL KIT

KTI-0002



| | | | | | | | |
|---|---|---|---|---|--|--|--|
|  |  |  |  |  |  |  |  |
| TD1-0002 | TD1-0003 | TD1-0005 | TD1-0006 | TD1-0007 | TD1-0008 | TD1-0009 | TD1-0010 |
| Marking Drill Ø 1.90mm | Pilot Drill External Irrigation Ø 2mm | Drill External Irrigation Ø 2.80mm | Drill External Irrigation Ø 3.20mm | Pilot Drill External Irrigation Ø 3.80mm | Drill External Irrigation Ø 4.50mm | Drill External Irrigation Ø 5.00mm | Drill External Irrigation Ø 5.50mm |
|  |  |  |  |  | | | |
| TD1-0019 | TD1-0020 | TK1-0012 | TK1-0013 | TK1-0015 | | | |
| Countersink Ø 3.75mm-Ø 4.2mm | Countersink Ø 5mm-Ø 6mm | Parallel Pin Ø 2mm, L 11mm Ø 3.20mm, L 11mm / S.S | Implant indicator L 14mm / S.S | Ratchet 84mm / S.S | | | |
|  |  |  |  | | | | |
| TK1-0002 | TK1-0004 | TK1-0010 | TK1-0019 | | | | |
| Key 1.28mm Long 23mm / S. S | Key 2.42mm Long 23mm / S. S | Motor key 2.42mm Long 28mm / S. S | 2.42 Retriever 15mm / S.S, Titanium | | | | |