



BY NOGA
MEDICAL

SOS Surgical kit.

MUST HAVE KIT FOR EVERY CLINIC

SCREW REMOVAL KIT

KT1-0009



TAG dental screw removal kit facilitates the removal of a broken screw in most cases with the possibility to renew the damaged thread.

Final user: All Doctors using implant internal hex (not only for TAG customers).



Screw Removal Burr Guide



1-72 UNF Tap



Screw Removal Burr

PLANNING KIT

KT1-0010



Easy to use planning kit allows you to choose the right abutment giving information on the required angulation and the height of the gingiva simultaneously.

Final user: All Doctors using implant internal hex. / Technicians. Recommended: TAG customers.



Mock Abutment Analog Straight



Mock Abutment Analog 15°



Mock Abutment Analog 15°

TECHNICIAN HOLDING KIT

KT1-0005



It's not always easy to reshape an abutment given its tiny size. This tool will allow you to easily maintain your abutment in order to give it the desired shape. It allows to work abutments with internal hexagon and conical connection abutments.

Final user: All Doctors using implant internal hex (not only for TAG customers).



1/4" Abutment Fixture
Long / 15 mm Stainless Steel



CC 1/4" Abutment Fixture
Long



Handle Set
100 mm Stainless Steel



1.28 Kit Screwdriver
140 mm / Stainless Steel

UNIVERSAL HOLDING KIT

KT1-0004



A simple tool that allows the surgeon to easily insert implants mainly in the anterior maxilla region.

Final user: All Doctors using implant internal hex (not only for TAG customers).



Handle Key 1.28mm/43mm
Stainless Steel



Extension for 1/4" Holder
50 mm Stainless Steel



Handle Key 2.42mm/43mm
Stainless Steel



CC Handle Key
2.42 mm



Handle Set 100mm
Stainless Steel